



Dermal Solutions Australia

Equipment Catalogue

**Australia's Leading Supplier - Professional Medical
& Beauty Equipment**

**CALL NOW FOR A FREE
NO OBLIGATION QUOTE
07 54791881**

**PROFESSIONAL MEDICAL
GRADE - BODIDERM
SKIN CARE PRODUCTS**

ABOUT US

WELCOME TO DERMAL SOLUTIONS AUSTRALIA – NO 1 IN PROFESSIONAL EQUIPMENT

Our Company strives to achieve only the best results from our Medical Grade Professional Equipment which should be your choice for effective in salon treatments. This catalogue provides you with the opportunity to view our range of products for Salon, Clinic and Medical Practices. We only provide equipment with the latest technology and highest standard in the world today. Our promise is excellent quality, value and proven results.

OUR MISSION

**Professional Clinical Results
Our Success Is Your Success**

OUR VALUES

Our professionally trained staff bring to you a lifetime of knowledge to provide you with the best training, support, marketing and medical grade equipment for your total success.

BEST PRICE – BEST SERVICE – BEST VALUE



Permanent Hair Removal Systems

There is a new demand for professional laser, IPL and now new SHR machines for better permanent hair removal. Clinic owners are now requiring state of the art Medical Grade equipment to provide this service to a wide range of customers. Skin types 1 – 6 can be treated without risk with the latest machines that Dermal Solutions offer. This will provide clients with fast, effective pain free hair removal systems that have better technology for smooth hair free treatments

SHR / SSR IPL Technology

New Innovation in IPL Systems – Laser Technology

SHR stands for Super Hair Removal, offering the latest technology to achieve long lasting hair removal. SHR uses photo selective thermolysis and combines it with

IN-Motion technology. The light in certain wave lengths and the heat emitted by SHR targets the pigment, the melanin, located in hair follicles. Flashes of light disable the targeted hair follicle, causing unwanted hair to fall out preventing future growth.

What makes SHR Stand out in the industry is its use of light pulse repetition rates at low energy levels, together with high power. Instead of simply blasting unwanted hair, it heats up the hair follicles very gradually without any discomfort, but maintains effectiveness so as not to damage the skins surface. The integrated cooling system causes no discomfort.

The SHR uses In- motion technology, a rolling technique that eliminates the common missed spots. The comprehensive coverage means smooth legs, arms and backs etc with the feeling of a hot stone massage.

Additional Features

The SSR module offers skin rejuvenation which reduces brown spots, capillaries, sun damage with no down time. The SSR technology repairs skin defects, treats wrinkles using low energy, short pulse, broadband light. Age spots, rosacea and imperfections can be eliminated.

Model: SHR SSR

Pulse Number: Single

Pulse Frequency: 3,4,5,6,10 hz

HR Spectrum: 690 – 950 nm

SR Spectrum: 560 – 950 nm

Energy Density: 1 – 36 j/cm² & up to 50 j/cm²

Hertz 1 – 10 Htz



Flexible Optional Extras

This machine can have RF added for skin tightening and firming as an added extra feature and handpiece. The machine can also come in SHR on one side and IPL on the other instead of SSR. This option includes filters 420nm, 530nm, 560nm, 640nm, 590nm, 755nm but these will need to be placed into the treatment head. Energy Density 1 – 50 Joules

SSR has a built in filter system for skin rejuvenation and pigmentation treatments so does not require manual filters. Energy density 1 – 36 Joules

SHR has a built in filter system for hair removal so does not require manual filters. Energy density 1 – 36 Joules

Faster & More Effective

The Dermalite SHR is a lot more powerful than normal IPL. It fires in Hertz with laser technology at a faster rate of delivery to the target areas. The built in treatment heads with own filters, means you do not have to insert small filters. That can be time consuming plus runs the risk of operator error if inserting an incorrect filter to a treatment head. Now you can relax and choose your skin type, working time and fluence more easily.

The treatment head can glide across the skin delivering faster more effective treatments.

Comfort

The crystal sapphire treatment head is cooled to minus 4 degrees during treatment so your clients will feel no pain and just a feeling of warmth. The machines system is water cooled, intercooled and air cooled providing the latest technology in cooling systems.

Reliability

German motors, Japanese Parts, American Xenon Lamps provide stable, effective and reliable results with the latest modular designs available.


Convenience

The Dermalite SHR is a very convenient table top version which can be put onto a beauty trolley and then moved from room to room. Most of our clients doing mobile work choose this machine for reliability, comfort and greater results.

Economical

The treatment heads are designed to last for up to 1 million flashes.

The multi system can treat all types of hair even red and blonde hair plus darker skin types from 1 – 6. It will treat vascular lesions, freckles, pigmentation, enlarged pores, and skin rejuvenation in a few short sessions.



DermaLITE Portable IPL

The Dermalite Portable IPL offers 3 different wavelengths for hair removal, 640, 690, 755nm for blonde, red and dark hair removal treatments. The large spot size of 15 x 50 offers operator ease in treatments with a very cold treatment head. The integrated cooling system is the best on the market today reaching -5 degrees for superior comfort and speed.

Applications

1. Hair Removal
2. Pigmentation Removal
3. Vascular Treatments – Spider Veins
4. Skin Rejuvenation
5. Acne Treatments

Specifications

Handpieces: 1 x IPL or 1 x SHR

Pulse Energy: 10 – 50 J/cm²

Pulse Sequence: 1 – 20 Pulses - adjustable

Pulse Delivery: 2 – 60 ms

Spot Size Options: 8 x 40, 12 x 40, 15 x 35 – (15 x 50 SHR Optional 690 – 950 NM)

Pulse Width: 2 – 20 ms

Filters: 420, 530, 560, 640 – 1200nm, 690 / 750 – 1200nm optional

Full Colour 10 x 4 cm Touch Screen

SHR Repetition Rate: 1- 10 HZ

SHR Shots: 1 – 60 seconds

SHR Energy: 1 – 36 jcm/2 or 1 – 50 jcm/2

Weight 38 kg

Operating Interface: New Pre Programmed Touch Screen (Pulses, Skin Type etc)

Options: IPL, SHR

Cooling Head: New Triple Cooling System For Treatment Head

Adjustable Cooling Head

Integrated Skin Cooling: -4 to 0 degrees

Cooling System: Water, Forced Air & Semi Conductor Cooling

American Xenon Lamp with 100,000 shots maximum, 300,000 shots for SHR

TGA, ISO, FDA, CE Medical Device Class IIA approvals

One Year Extended Warranty Parts Replacement

Client & Operator Goggles

Manual & Large Carry Case 0.4" intelligent touch screen interface system, easier to operate with standard and professional mode



- Filters coated by advanced film plating to make it more effective
- Self-inspection system
- Light protective handpiece – will not shoot if filter inserted incorrectly
- Every treatment method has high performance from extensive clinical studies
- Accurate energy output, to meet ideal operating environments
- Advanced integrated cooling system with water & air & semiconductor
- Modular design -easier to check, upgrade and maintain.

Cavitation RF Systems

CAVITATION

Is a Non-Invasive Treatment system for the safe reduction of adipose tissue and cellulite. It works using a multi-frequency mode that permits a careful, selective release of ultrasound waves combined with a process called "Cavitation" which selectively targets fat cells in specific problem areas to safely disrupt the membranes of unwanted fat cells without affecting the surrounding tissue such as nerves or blood cells.

WHAT IS CAVITATION?

Cavitation is a biological phenomenon that consists of the formation of pockets of steam (Bubbles) within the interstitial (in between cells) areas which implode by the action of a pulsed frequency ultrasound. This creates shock waves that selectively damaging the membranes of the adipocytes (fat cells), spilling the liquid fat into the blood stream.

WHAT HAPPENS DURING TREATMENT?

Treatments are extremely relaxing, and pain-free and do not last any longer than 20-40 minutes per visit. The client relaxes on the treatment bed while a hand held flat probe device is glided over the marked areas delivering the multi-frequency Cavitation energy. Following the treatment, the client can immediately resume their daily activities.

Treatments require no anesthetics, no downtime, a maximum of 25 minutes (Excluding consultation time) and gives you a quick and a worthwhile return on your investment.

CELLULITE is a change in the fat stored under the skin, appearing as dimples or ripples in the skin. It affects most women. Studies have shown that 80 % of women have cellulite, even if they are thin.

The dimples can go through stages, actually getting worse until the skin sags and loses collagen and elastin. Our Bodilite machines will blast the fat cells, shrink the collagen and therefore reduce cellulite and tighten the skin.

Focal heating of the fat cells increases the rate at which the liquid fat is released from the tissue. It will not affect the surrounding tissue but will increase the secretion of fat.

RF Systems

RF Skin Tightening is a new technique for skin lifting and contouring with no surgery. Amazing results are seen after one or two treatments and can last for up to 5 years. *RF treatments* are suited for the people who may have experienced some skin sagging and lack of firmness resulting in looser and wrinkled appearance.

RF Skin Tightening works to generate new collagen for up to 6 months after the initial treatment and to tighten the deeper dermis resulting in an improvement in skin texture, tone and contour. *RF skin treatment* is perfectly suited to treat the sagging of the jawline, face, eyes, tummy arms and legs. There is no downtime and treatments are painfree.

How does it work?

Rf skin tightening is a radio frequency sound wave that is safe and non-surgical that gently heats the collagen fibres immediately tightening them. New collagen is stimulated by the fibroblasts in the deeper layer of the skin over the next few months. The radio-frequency also improves blood flow and lymphatic drainage thereby encouraging a better oxygen supply and toxin removal to further reduce the effects of ageing.

How does the treatment feel?

The gentle warmth experienced during treatment is similar to a hot-stone massage – very soothing and pain free. There is an immediate improvement of the lax or loose skin and the improvement continues for up to 6 months after the treatment.

Do I need more than one treatment?

RF treatment is a gentle non-surgical procedure and as such your technician will gently treat your problem areas. It is advisable to have a few treatments over a 6-8 week period for the ultimate results to allow the deeper layers to rejuvenate and regenerate.

Versatile

This machine comes in the multi platform system with RF Skin Tightening plus the added bonus of a vacuum liposuction handpiece on the portable model for lymphatic drainage. RF and Cavitation can be added to the multiplatform systems.



Laser ND Yag Tattoo Removal

The laser device contains two treatment heads: 1064nm & 532nm. The 1064nm head is used for blue or black pigmentation and the 532nm head is for the red, yellow and brown pigmentation.

The machine will treat tattoos, permanent makeup, freckles and yellow / brown / coffee spots, birthmarks and age spots.

There is a water cooling system in the laser device for cooling down the laser light for client comfort.

The soft ware of the laser has two main functions, to allow the operator to select the best individual parameters for each patient.

Table Top Version

- Laser Type: Q switched Nd YAG Laser
- Screen: 12.4 inch full colour touch screen
- Laser wavelength: 1064nm & 532nm
- Pulse energy: 0 - 900mj & 0- 300mj
- Pulse width: <12ns
- Repeat frequency: 1、2、3、4、5、6HZ (0.5 – 3 Hz optional)
- Light conducting method: conducted directly from the laser outlet
- The light spot diameter of the laser output 1 – 6 mm adjustable
- Environment temperature: 5°C~40°C
- Relative humidity: ≤70%
- Weight 25 kg
- Cooling system: Installed circulating water cooling+ air cooling.



Features

Latest technology, energy output is higher, no need to calibrate manually, more reliable and longer life time

Accurate treatments that will not damage normal tissue

No downtime

Significant performance without side effects

Accurate energy output with a safe operating environment

Advanced integrated cooling system with water, air & semiconductor

Module design - easier to check - upgrade and maintain


Large Stand Up Version

- Laser Type: Q switched Nd YAG Laser
- Screen 12.4 inch full colour touch screen
- Laser wavelength: 1064nm & 532nm
- Pulse energy: 0 - 1200mj & 0 - 400mj
- Pulse width: <10ns
- Repeat frequency: 1、2、3、4、5、6HZ (0.5 – 3 Hz optional)
- Light conducting method: conducted directly from the laser outlet
- The light spot diameter of the laser output 1 - 7mm adjustable
- Environment temperature: 5°C~40°C
- Relative humidity: $\leq 70\%$
- Cooling system: Installed circulating water cooling+ air cooling.
- Weight 48 kg
- Unique design – Quality Handpieces with removable tips



Laser Features

Fragments Ink – Will Not Damage Surrounding Tissue



Unique Modern Designed Handpieces -
Operator Speed & Comfort



Multi Function Machines

Powerful multi functioning devices allow a business to have more treatments in just one machine. Add IPL, Ultrasonic Fat Cavitation, SHR & RF Skin Tightening to your machine of choice. Ask today about the benefits of all treatments and expand your client service menu.

More Treatments, Happier Clients, Proven Treatment Options Now Available

Applications

1. Hair Removal
2. Pigmentation Removal
3. Vascular Treatments – Spider Veins
4. Skin Rejuvenation
5. Acne Treatments
6. Ultrasonic Fat Cavitation 40,000 Hertz and 1 Mega Hertz Heads
7. RF Skin Tightening
8. SHR – Super Hair Removal 690 – 950 nm laser technology
9. ND Yag Laser Tattoo Removal
10. IPL – Hair Removal & Skin Rejuvenation

Pain Free

The triple cooling system plus cold treatment heads to minus 4 degrees offer you total client comfort and confidence.

Versatile

With so many options available, you can tailor make your machine to suit your individual clinic needs. Patented Technology provides you with stable, effective treatment options.

Reliable

The systems offer the latest up to date technology providing you with a multi platform system with total reliability in performance and parameters.

Convenience

The Multi Platform machine is on wheels and can be moved easily from room to room.

HYDRODERMABRASION MICRO EQUIPMENT

Hydra Peel 4 in 1 The machine has superior strength in application combining this with a blasting of water onto the skin with the hydrodermabrasion wand and tips. The vacuum pressure and water pressure can be adjusted individually to suit individual client needs. The diamond dermabrasion wand has 12 different tips in sizes and abrasiveness for thickened, sun-damaged skins and scarring. Water is stored in an airtight sealed container and it has a disposable sealed water container for used water. An added feature is an alarm to let the operator know when the water bottle is near empty.

The oxygen spray feature can be used post treatment to calm and soothe inflamed skin plus the operator can add serums to the bowl on the wand to infuse products.



Hydrodermabrasion utilizes the natural healing powers of water and oxygen to effortlessly exfoliate the skin without the use of hard crystals or abrasive textured wands, generating deeply hydrated and healthier looking skin. Dermal Solutions has created a crystal free dermabrasion applicator, a hand held nozzle that works by utilizing a high velocity jet of water combined with oxygen. This completely natural combination of micro droplets of water is dispersed from the nozzle at high speed, striking against the skin and thereby removing any dead skin cells with little effort and absolutely no pain to the client, deeply exfoliating and rejuvenating the skin.

Hydrodermabrasion is an amazing step forward in skincare technology with highly effective results seen after just one treatment

The Diamond Tip Microdermabrasion comes with seven diamond tips and two stainless steel handpieces, filters plus more. The tips are stainless steel for easy cleaning and can be used on the face and neck or anywhere on the body to remove dead skin cells. Each tip is a different size including a smoothing tip to vacuum up debris on completion. Ask us today about the best machine to suit your clinic's needs. The new oxygen spray with attachments to infuse products adds the most amazing finish to your treatments.

Technical Data

Electricity voltage: AC90V-220v

Working frequency: 50Hz/60HZ Out put: 6 kpa

Power: 250w

Packing size: 62*53*34cm

Gross weight: 25KG

Accessories list accessories

Main body	1pcs
Clean water bottle	1pcs
Waste water bottle	1pcs
Hydro handle	1pcs
Hydro facial tips	8pcs
Micro-dermabrasion handle	2pcs
Micro tips	9pcs
Oxygen spray handle	1pcs
Black rings	3pcs

SKIN NEEDLING

BODIDERM has a new titanium Derma Roller which is a revolutionary hand-held device. It has the ability to naturally stimulate collagen and elastin in the skin to promote healthier skin cells in the body. The fine titanium medical-grade needles on the dermaroller are painless. They create micro channels in your skin which accelerate the body's natural healing process. This creates new healthier tissue to replace the imperfections in the skin. As it renews your skin cells regularly it will treat common problems such as wrinkles, skin hyper pigmentation, acne scars, uneven skin and stretch marks.



What can the Derma Roller System treat?

- Anti Ageing
- Anti Wrinkle
- Hair Loss Treatment / Hair Restoration
- Hyper Pigmentation Treatment
- Scar Removal - including acne scar removal
- Skin Smoothing
- Stretch Mark Removal

When we try to *resurface* we target the epidermis (top layer of skin). When we need to remodel the skin for wrinkle treatments, scar treatments we are targeting treatments within the mid-dermis. The dermarollers' fine needles create small micro-channels in the skin leaving the epidermal barrier fully intact, treating only areas in the mid-dermis.

Using the Derma Roller System with Hydrasilk serum allows Hyaluronic acid and healing minerals to penetrate through the top layer of the epidermis through the micro-channels and repair the skin from the inside out without damaging the surface of the skin. This eliminates the risks and negative side effects on the skin commonly encountered other invasive skin peels.


While your pores are opened the dermaroller stimulates collagen and elastin production. The production of collagen and elastin helps to thicken the skin, smooth away fine lines and wrinkles and rid the skin of hyper pigmentation. When the dermaroller is applied simultaneously with the application of Hydrasilk your skins absorption of vitamins, enzymes, lipids and other nourishing ingredients are greatly enhanced. In approximately one to two hours after the application of the dermaroller along with Hydrasilk these micro-channels fully close through the body's natural healing process, providing just enough time for new collagen structures to form within the skin's lower layers.

The Derma Roller System is effective in acne scar removal as well as treating breast augmentation scars, burn scars, C-section scarring, chicken pox scars, ice pick scars, large pores, pitting of the skin, surgical scars, tummy tuck scars plus more.

The Derma Roller System can be used to remove stretch marks such as those developed during pregnancy, adolescent growth spurts, rapid muscle growth, rapid muscle growth spurts and weight loss/gain.

The Rollers come in sizes:

0.20mm, 0.25mm, 0.50mm, 0.75mm, 1.00mm, 1.50mm, 2.00mm, 2.5mm.



DERMAL SOLUTIONS PRICE LIST – PLEASE ADD 10 % GST

IPL SYSTEMS:

PORTABLE IPL MACHINE	\$22,000.00
PORTABLE IPL WITH ELIGHT / RF	\$25,500.00
STAND UP ULTRA IPL – 2 HEADS	\$43,000.00
STAND UP IPL ULTRA WITH ELIGHT / RF	\$47,000.00
NEW SHR / SSR/ RF (Laser Technology)	\$36,000.00
NEW SHR / SSR ROLLING (Laser Technology)	\$33,000.00

CAVITATION SYSTEMS:

ULTRA CAVITATION WITH RF / VACUUM	\$21,500.00
ALL IN ONE CAVITATION / IPL / RF	\$49,000.00

MICRODERMABRASION:

DIAMOND TIP	\$3,500.00
HYDRODERM / DIAMOND TIP/ OXYGEN	\$8,500.00

LASER EQUIPMENT:

ND YAG QSWITCHED TATTOO REMOVAL	\$30,000.00
ND YAG QSWITCHED TATTOO / VEINS / SKIN	\$33,000.00
ND YAG PICOSURE LASER TATTOO REMOVAL	\$75,000.00
ND YAG 808 NM HAIR REMOVAL	\$65,000.00
FRACTIONAL CO2 LASER SKIN REJUVE	\$68,000.00
SKIN TAG / SPIDER VEIN REMOVAL	\$12,000.00
PDT LED MACHINE	\$ 8,500.00

CONSUMABLE EQUIPMENT:

CLIENT GOGGLES	\$ 85.00
OPERATOR GOGGLES IPL	\$ 220.00
OPERATOR GOGGLES LASER	\$ 380.00
IPL FILTERS	\$ 220.00
IPL HANDPIECE (CHOICE OF SPOT SIZE)	\$2,500.00
SHR / SSR HANDPIECES	\$2,800.00
ND YAG LASER HANDPIECES	\$2,500.00
CAVITATION & RF HANDPIECES	\$ 800.00
WATER FILTER EXTERNAL	\$ 75.00
LASER TIPS 532NM & 1064NM	\$ 500.00
DERMAROLLER	\$ 75.00



Dermal Solutions Australia

DERMAL SOLUTIONS AUSTRALIA
SHOP 3 / 31 – 33 PLAZA PARADE,
MAROOCHYDORE 4558 QLD PH: 0754791881
WEBSITE: www.dermalsolutions.com.au