Transducer Guide LOGIQ™ P9 and LOGIQ P7





	Description	Application	FOV	Bandwidth	Biopsy Guide	System
	Convex					
C1-5-R H40462L		Abdomen, OB/GYN, Urology	70°	1 – 6	Multi-Angle, disposable with a reusable bracket (H40432LE)	LOGIQ P9
4C-R H4000S		Abdomen, OB/GYN, Urology	55°	1 – 5	Multi-Angle, disposable with a reusable bracket (E8385NA)	LOGIQ P7
	Micro-convex					
8C-R H40402L		Neonatal, Pediatrics	132°	3 - 11	No	LOGIQ P9 LOGIQ P7
E8C-R H40402L	-	OB/GYN, Urology, Endocavity	123°	3 - 11	Single-Angle, disposable with a disposable bracket (E8385MJ, E8333JB), Single-Angle, Reusable bracket (H40412LN)	LOGIQ P9 LOGIQ P7
	Linear					
9L-R H40442L		Vascular, Small parts, Pediatrics, Abdomen	44 mm	2 – 8	Multi-Angle, disposable with a reusable bracket (H4906BK)	LOGIQ P9
12L-R H40402L		Small parts, Vascular	37 mm	3 – 12	Multi-Angle, disposable with a reusable bracket (H40432LC)	LOGIQ P9
L6-12-R H48062A		Small parts, Vascular, Pediatrics	37 mm	4 - 13	Multi-Angle, disposable with a reusable bracket (H40432LC)	LOGIQ P7





© 2014 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

 $\ensuremath{\mathsf{GE}}, \ensuremath{\mathsf{GE}}$ monogram, and LOGIQ are trademarks of General Electric Company or one of its subsidiaries.

GE Medical Systems Ultrasound & Primary Care Diagnostics, LLC, a General Electric company, doing business as GE Healthcare.

GE Healthcare 9900 Innovation Drive Wauwatosa, WI 53226 U.S.A. www.gehealthcare.com